

IMPORTANT INFORMATION CONCERNING THE MAGNUM FOOT

Please read thoroughly before fabricating a prosthesis using the Magnum Foot.

Warranty

The warranty for the Magnum Foot is one year (allowing three weeks from the date of purchase for installation), provided that the foot is selected according to the following criteria:

Note: The weights listed below are adjusted body weights. Adjusted body weight is defined as the weight of the amputee plus any loads carried by the amputee.

	Level 1	Level 2	Level 3	Level 4
Magnum Foot	350 lbs 159 kg	350 lbs 159 kg	350 lbs 159 kg	300 lbs 136 kg

Level 1: Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Level 2: Has the ability to traverse low level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the community ambulator.

Level 3: Has the ability or potential for ambulations with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Level 4: Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Use of a Magnum Foot for amputees who do not meet the above criteria or who engage in extremely high and abusive activities is against Ohio Willow Wood's recommendations and will void the one-year warranty. "Extremely high and abusive activities" are defined as activities such as skydiving, karate, and judo; activities that could result in injury to an individual's natural limbs; and activities that expose the prosthesis to corrosives such as salt water.

Modifying the Foot

The urethane foam shell may be lightly sanded with little resulting loss of strength. Do not expose the internal components of the foot, as this could reduce the structural integrity of the foot.

Exposure to Water

The Magnum Foot is essentially water resistant, provided no open cracks have developed and no area of the foot has been sanded. Complete submersion in water is not recommended.



OHIO WILLOW WOOD®

free the body...free the spirit®

15441 Scioto Darby Road

Mt. Sterling, OH 43143

phone 740.869.3377 / 800.848.4930

fax 740.869.4374 www.owwco.com



Ohio Willow Wood Company B.V.
Keizersgracht 62/64
1015 CS Amsterdam
The Netherlands



PN-1477-D 9 NOV 2010

U.S. Patent #5,728,171

Patient Advisory Warning

The attached Patient Advisory Warning enables you, the prosthetist, to effectively notify your patients of the limitations of the components in their prosthesis, and of the need to monitor their weight and activity levels. Please review the Patient Advisory Warning with the patient upon delivery of a prosthesis with a Magnum Foot. The patient and the prosthetist should then sign the Patient Advisory Warning to acknowledge that it has been reviewed and understood by both parties. Give one signed copy to the patient and place one copy in the patient's file.

If a patient's weight or activity level increases after receiving a prosthesis with a Magnum Foot, the patient should immediately contact the prosthetist to determine whether replacement components are necessary. If a patient continues to use a prosthesis with a Magnum Foot after experiencing an increase in weight and/or activity level, the foot could fail with the possibility of serious injury to the patient.

To ensure that the correct components are selected for each patient, the prosthetist should weigh the patient on scales in the prosthetist's office. Do not rely on the patient's estimate of his/her own weight. Instruct the patient to monitor his/her weight weekly to ensure that it remains in a range appropriate for the prosthetic components being used.
