


INSTALLATION INSTRUCTIONS FOR THE MAGNUM TITANIUM SLIDE ADAPTERS

Please read instructions thoroughly before installing the Magnum Titanium Slide Adapter with 4-Hole Rotatable Pyramid (MHD-224150) or the Magnum Titanium Slide Adapter with 4-Hole Pyramid (MHD-224125).

Note: The pyramid dome on the Magnum Titanium Slide Adapter with 4-Hole Rotatable Pyramid is held in place by a magnet with an intensity of 28 gauss at the magnet surface, and 0 gauss at a distance of 4" from the surface. This intensity is comparable to that of a refrigerator magnet. Please share this information with the amputee, in the event that the amputee is wearing a pacemaker or any other device that could be affected by magnetism.

Recommended Installation Procedures

1. Using the four M6 X 16 screws included with the Slide Adapter, attach the Slide Adapter to the socket as shown. (Alternately, use the M6 X 35 screws that are also included with the Slide Adapter to attach the Slide Adapter to an Alpha Lock or Alpha Locking Lanyard.) Tighten the screws to 9 ft-lbs (12 Nm).



2. Connect a pyramid receiver component of your choice to the Slide Adapter. Tighten the setscrews using the torque and threadlocker specifications provided by the manufacturer of the pyramid receiver.

Note: the purpose of the 5mm setscrew is only to keep the pyramid from sliding when angular adjustments are made. If the 5mm setscrew is tight before attaching the receiver, it will prevent the dome from locking tightly against the housing and titanium pyramid.

3. To make a slide adjustment:
 - a. Loosen the 5mm setscrew on the side of the Slide Adapter.
 - b. Loosen two adjacent 8mm setscrews on the pyramid receptacle.
 - c. Slide the pyramid to the desired position.
 - d. Before tightening the 5mm setscrew, tighten the 8mm setscrews using the torque and threadlocker specifications provided by the manufacturer of the pyramid receptacle.
 - e. Hand-tighten the 5mm setscrew on the side of the Slide Adapter.
4. To make a rotation adjustment (applicable only if using MHD-224150):
 - a. Loosen two adjacent 8mm setscrews on the pyramid receptacle.
 - b. Rotate the pyramid to the desired position.
 - c. Tighten the 8mm setscrews using the torque and threadlocker specifications provided by the manufacturer of the pyramid receptacle.



OHIO WILLOW WOOD®

free the body...free the spirit®

15441 Scioto Darby Road

Mt. Sterling, OH 43143

phone 740.869.3377 / 800.848.4930

fax 740.869.4374 www.owwco.com

PN-1653-D 26 MAR 2007

Warranty

The warranty for the Magnum Titanium Slide Adapters is two years, provided that the adapters are selected according to the following criteria:

Note: The weights listed below are adjusted body weights. Adjusted body weight is defined as the weight of the amputee plus any loads normally or routinely carried by the amputee.

	Level 1	Level 2	Level 3	Level 4
Magnum Titanium Slide Adapter with 4-Hole Pyramid	350 lbs 159 kg	350 lbs 159 kg	350 lbs 159 kg	300 lbs 136 kg
Magnum Titanium Slide Adapter with 4-Hole Rotatable Pyramid	350 lbs 159 kg	350 lbs 159 kg	350 lbs 159 kg	300 lbs 136 kg

Level 1: Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Level 2: Has the ability to traverse low level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the community ambulator.

Level 3: Has the ability or potential for ambulations with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Level 4: Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Use of a Magnum Titanium Slide Adapter for amputees who do not meet the above criteria or who engage in extremely high and abusive activities is against Ohio Willow Wood's recommendations and will void the two-year warranty. "Extremely high and abusive activities" are defined as activities such as skydiving, karate, and judo; activities that could result in injury to an individual's natural limbs; and activities that expose the prosthesis to corrosives such as salt water.

Warranty Disclaimer

Ohio Willow Wood warrants that each product manufactured will, at the time of delivery, be of workmanlike quality and substantially free of defects. OHIO WILLOW WOOD MAKES NO OTHER WARRANTY, IMPLIED, OR EXPRESSED, AND MAKES NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. This warranty shall terminate immediately upon an action to combine our products with other materials or in any manner to change the nature of our products. The sole remedy is replacement of the products or credit for the products. Ohio Willow Wood's liability shall not exceed the purchase price of the product. Ohio Willow Wood shall not be liable for any indirect, incidental, or consequential damage.

Ohio Willow Wood Retention of Rights

Ohio Willow Wood retains all intellectual property rights reflected or incorporated in its physical products, regardless of the transfer of the physical products to another party or parties.
