LimbLogic® M Diagnostic Kit Instructions


WHAT'S IN THE BOX

Diagnostic Plate	Tube
Proximal Gasket	Vacuum Gauge
M6 X 60 mm Socket Head Capscrews (4)	Instructions
Diagnostic Disc	

The Diagnostic Kit can be used to evaluate the vacuum level of the system to diagnose the source of a leak.

1. Remove the pump from the socket.

2. Place a gasket on the Diagnostic Plate.

3. Place the Diagnostic Plate with the gasket on top of the pump inlet and gasket.

*two gaskets:
one on top
of Diagnostic
Plate;
one between
pump and
Diagnostic Plate*

4. Bolt the pump onto the socket with the M6 x 60 mm bolts.

5. Install the gauge.


Ask the amputee to activate the pump by shifting their weight on and off of the prosthesis. Observe the vacuum level with the gauge.

- If the vacuum level is adequate, remove the Diagnostic Plate and reinstall the pump with the M6 x 50 mm bolts.

• If a leak is detected:

a. Remove the pump. Install the Diagnostic Disc as shown.

b. Bolt the pump onto the socket with the M6 x 60 mm bolts.

c. Install the gauge.

Ask the amputee to activate the pump by shifting their weight on and off of the prosthesis. Observe the vacuum level with the gauge.

- An adequate vacuum level indicates that the leak is in the sleeve or the socket.
- If a pump leak is detected, call WillowWood Customer Care at 800.848.4930.


6. Remove the Diagnostic Plate and Diagnostic Disc before allowing the patient to leave with the prosthesis.

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Ohio Willow Wood Company B.V.