

REFURBISHED LimbLogic® VS

- All repairs performed by WillowWood
- Refurbishment fees not discounted
- Units thoroughly tested before shipping to customers

Option 1: Basic Refurbishment

Indications	System is not working consistently or will not operate.
Includes	Filters, tubing and check valves will be replaced. These are a combination of minor components within the unit that if loose or clogged can impair or prevent correct operation of the LimbLogic VS unit. The system is fully tested before being returned to the customer.
Warranty	3 months from invoice date
Lead Time	5 business days from date of refurbishment order

Option 2: Advanced Refurbishment

Indications	System is not working consistently or will not operate.
Includes	Filters, tubing, check valves, and battery will be replaced. In addition, the internal electronics and pump of the LimbLogic VS unit are replaced. The system is fully tested before being returned to the customer.
Warranty	6 months from invoice date
Lead Time	5 business days from date of refurbishment order