


# WillowWood® Express Liner Sizing Guide

Please refer to the illustrations for recommendations for fitting WillowWood® Express Liners. Shown below are limb circumferences that the various liners will fit, measured at 4 cm (1.5 in) and 30 cm (12 in) from the distal end.

If converting from Alpha® Classic or Alpha Hybrid Liners with a uniform profile, use the two measurements to determine the appropriate size in WillowWood Express Liners.



	Medium	Medium Plus	Large	Large Plus
<b>Proximal Measurement</b>	23 - 36 cm (9 - 14 in)	33 - 50 cm (13 - 20 in)	36 - 55 cm (14 - 21.5 in)	40 - 60 cm (15.5 - 23.5 in)
<b>Distal Measurement</b>	18 - 26 cm (7 - 10 in)	20 - 28 cm (8 - 11 in)	25 - 33 cm (10 - 13 in)	28 - 40 cm (11 - 15.5 in)

A diagram showing four trapezoidal liner profiles. From left to right, they are labeled 'Medium', 'Medium Plus', 'Large', and 'Large Plus'. Each profile has a dashed line indicating the proximal measurement (30 cm from distal) and another dashed line indicating the distal measurement (4 cm from distal). The profiles are shaded to show depth.

	Medium	Medium Plus	Large	Large Plus
<b>Proximal Measurement</b> 30 cm (12 in) from distal	M	M+	L	L+
<b>Distal Measurement</b> 4 cm (1.5 in) from distal				