

WillowWood® Custom Solutions Order Form

WillowWood One Transfemoral System with LimbLogic® Definitive Socket

GENERAL

Date _____ P.O. number _____
Prosthetist name _____ Purchased From Cascade SPS
Email address _____ WillowWood Other _____
Facility _____ Account number _____
Ship to _____
Shipping address _____
City _____ State _____ ZIP Code _____
Phone _____ Preferred Contact Method _____
Scheduled Patient Fit Date _____

PATIENT INFORMATION

Patient ID: _____
Amputation side Right Left
Special instructions: _____

Items sent with Socket: _____

Please place LimbLogic Pump
Serial Number label here

*If a replacement socket, please record
LimbLogic pump serial number*

SOCKET INFORMATION (PLEASE CHECK)

Includes fabrication of positive model, trim, and finish.

Is this a replacement WillowWood One socket?
 Yes No

Pump Placement: A P M L

Valve Placement: A P M L

** Indicate Distal Adapter for use in Definitive Socket:

** Distal mount applications require the included LimbLogic distal adapter be used for definitive socket. For side mount applications, the distal adapter used in diagnostic socket will be used in definitive socket unless directed otherwise.*

NOTE: Changing Distal Adapter will incur additional charges.

**Be sure to indicate the alignment lines
(sagittal and coronal planes) on the socket.**

Lower Socket to Adapter? Yes No

Back Alignment Button: Exposed Covered

Color of Lamination: _____

**** The diagnostic socket should be the final modified version, in dynamic alignment, and with a distal adapter attached. WillowWood does not have the ability to modify sockets upon receipt.**

DO NOT SEND ANY OF THE FOLLOWING TO WILLOWWOOD: (Store in original box with Patient ID)

Knee Assembly	One Seal	Pump	Charger
Alpha® Liner	One Gel Sock	Fob	

SHIPPING METHOD (CHECK ONE)

Methods other than Ground are subject to additional charges.

FedEx Ground FedEx Priority Overnight
 FedEx Express 2-day A.M. Other _____

*Please retain a copy for your records.
Download a writable PDF at willowwoodco.com
in the Custom Fabrication section.*

**For order processing, please ship completed
order form and items to the nearest location:**

Utah Facility

1760 West Redwood Depot Lane, Suite 6
Salt Lake City, UT 84104 USA

Ohio Facility

15441 Scioto Darby Road
Mt. Sterling, OH 43143 USA

WillowWood®

Phone: 800-848-4930; Fax: 888-878-4858
www.willowwoodco.com