

OMEGA[®] Software

Version 1.4.4.4 Release Notes

OMEGA Software, version 1.4.4.3 is an update released in **August 2017** which includes enhanced features and fixes for reported issues.

For detailed information on using any of OMEGA's features and capabilities, select **Help** from the Main Menu Bar of the software. **Help** provides a link to the User's Guides as well as links for contacting WillowWood for advice or help.

Additional OMEGA product information can be found in the [OMEGA section of willowwoodco.com](http://www.willowwoodco.com).

1. ADDITIONS An alert will display for Windows 10 users who select software options that require scanning. The alert informs users that scanning is not available.

A **View Release Notes** button has been added to the **Available Update** window. The **View Release Notes** button allows users the opportunity to view the release notes prior to installing software updates.

2. REMOVALS

The **Ramus Containment Modification** tool has been removed from the software. However, this tool remains available in user-created Design Stack templates and in shape files that were modified and saved using this function.

3. CHANGES

Danmar order forms have been updated within the software.

4. FIXES

Bug fixes have been made to prevent software errors and to improve the use of the **TF by Measurements** feature.

WillowWood[®]

The Ohio Willow Wood Company
15441 Scioto Darby Road
Mt. Sterling, OH 43143
phone 740.869.3377 / 800.848.4930
fax 740.869.4374 www.willowwoodco.com



Ohio Willow Wood Company B.V.
Keizersgracht 62/64
1015 CS Amsterdam
The Netherlands
Patent www.willowwoodco.com/education-and-resources/patents

