

# Alpha DESIGN® Liner Ordering Instructions

Want to order a DESIGN Liner? Here's what you need to know.

## 1. SELECT ONE OF THE FOLLOWING SHAPE CAPTURE OPTIONS

- Use **OMEGA® Tracer®**.
- Use the [Alpha DESIGN App](#) with an iPad and a Structure Sensor.
- Submit a cast or scan with a completed [online order form](#) from our website.
- Submit a cast or scan via FusiformCAST. (Click [here](#) to request your invitation!)

## 2. CASTING (IF APPLICABLE)

Please follow the [BK casting instructions](#) or [AK casting instructions](#) on the website.

### Please Note:


- For bilateral amputees, please use a separate order form for each residual limb to achieve the best socket/liner accuracy.
- When casting transtibial patients, please cast at least 4" to 5" above the knee.
- Please write the patient ID on the cast.

Ship casts to: WillowWood ATTN: DESIGN Dept.  
15441 Scioto Darby Road  
Mt. Sterling, OH 43143 USA

## 3. ORDERING

Make sure to have the following information ready to enter into the software or the online order form:

- Measurements shown at right
  - Purchase Order number
  - WillowWood account number
  - All patient and liner information
  - Ship-to address
  - .stl or .aop file (if applicable)
- Note: Shape files can be uploaded to the online order form, or emailed to [aopfiles@owwco.com](mailto:aopfiles@owwco.com)**



After submitting the completed online order form, you will receive a summary of your order via email. **Please print the summary and include it in the shipping box if you are sending a cast.**

**Failure to complete and submit an online order form could delay your order.**

# WillowWood®

The Ohio Willow Wood Company  
15441 Scioto Darby Road  
Mt. Sterling, OH 43143  
phone 740.869.3377 / 800.848.4930  
fax 740.869.4374 [www.willowwoodco.com](http://www.willowwoodco.com)



Ohio Willow Wood Company B.V.  
Keizersgracht 62/64  
1015 CS Amsterdam  
The Netherlands  
Patent [www.willowwoodco.com/education-and-resources/patents](http://www.willowwoodco.com/education-and-resources/patents)

